
P

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en - Refer to Instructions for Use
fr - Consulter la notice

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INTERCEPT PROCESSING SET for PLASMA

For use with INTERCEPT Illuminator

Each set is wrapped in a tamper-evident package and includes one 15mL 6 mM amotosalen hydrochloride solution container (Formula : Amotosalen HCl 203 mg + Natr. chlorid. 924mg - Aqua ad iniect. ad 100 mL), one illumination container, one compound adsorption device (CAD), three INTERCEPT Plasma storage containers. The set is sterilized by a combination of steam and radiation.

Indications and Usage

This set is used with an INTERCEPT Illuminator to inactivate a broad spectrum of viruses, bacteria, parasites as well as contaminating donor leukocytes in plasma. INTERCEPT plasma is indicated for support of patients requiring plasma transfusions, according to clinical practice guidelines. INTERCEPT plasma may be stored from the time of collection for 12 months between -18°C and -25°C or for 24 months below -25°C, in compliance with applicable procedures and regulations.

Contraindications

Use of INTERCEPT Plasma is contraindicated in patients with a history of allergic response to amotosalen or psoralens.

Precautions

Do not use if : tamper-evident package has been opened; signs of deterioration are visible; fluid path closures are loose or not intact; cannulae are broken or there is no fluid in amotosalen solution container. Do not store above 25°C. Do not vent. Do not freeze. Protect the pack and tubing from sharp objects.

Unused sets in open aluminium foil may be kept 20 days at room temperature by folding and securing open end of aluminium foil. Units removed from the aluminium foil must be used within 24 hours.

Keep set in light-protective package until time of use. Protect from direct sunlight and strong UVA light source.

Set is single use only. Do not reuse. Do not resterilize. This product is not designed for reuse. Misuse can result in adverse reactions, including severe illness and possibly death.

All the following conditions must be met for pathogen inactivation:

- Plasma volume and red blood cell (RBC) content must be within the range specified in **Table 1**.
- Plasma mixed with amotosalen must be exposed to UVA light dose from INTERCEPT Illuminator. No other source of UVA light may be used.
- Plasma must be passed through the CAD by gravity flow process after illumination.
- Apheresis plasma must be frozen within 18 hours after collection. Pooled whole blood derived plasma must be frozen within 19 hours after collection of the oldest unit in the pool.

This process is designed to be a closed system. Treatment with INTERCEPT Blood System does not replace applicable standards for processing in open and closed systems. If there is a leak in the set during processing, plasma product must be discarded.

Warnings : Amotosalen in contact with skin may result in photosensitisation in the presence of ultraviolet light. If skin exposure occurs, flush exposed skin copiously with water. Sterile connecting device (SCD) should be used according to manufacturer's instructions for use.

Notes to Physicians

While laboratory studies of amotosalen processing with UVA light have shown a reduction in levels of certain viruses, bacteria and parasites; there is no pathogen inactivation process that has been shown to eliminate all pathogens.

INTERCEPT plasma components should not be prescribed to newborns treated with phototherapy devices that emit a peak energy wavelength less than 425 nm, and/or whose emission bandwidth limit is less than 375 nm, due to the risk of erythema resulting from potential interaction between ultraviolet light (below 400 nm) and residual amotosalen. A maximum daily intake limit of 1300µg amotosalen is recommended.

Instructions for Use

Materials Needed : One (1) INTERCEPT Processing Set for Plasma.

Equipment Needed : INTERCEPT Illuminator, Sterile Connecting Device (SCD), Tube Sealer.

Temperatures should be controlled to meet applicable regulations for plasma processing.

Used and unused INTERCEPT sets should be discarded like any used blood containers, as biohazardous waste.

Process Steps

A-Plasma for Processing with the INTERCEPT Blood System

Single donor or pooled plasma products within the following ranges have been shown to be acceptable for use with this processing set.

Table 1

Volume	RBC content
385-650 mL	<4 × 10 ⁹ RBC/mL

B- Amotosalen Addition

1. Remove set from package.
2. Weld tubing from plasma container to amotosalen container tubing using SCD.
3. Disassemble set from organizer and remove rubber band.

4. If two plasma units will be produced by the INTERCEPT process, heat seal and remove one storage container.
5. Label set containers using appropriate donation identification . See warning in section C. While labeling storage containers, separate them to ensure they do not adhere to one another.
6. Hang plasma container, ensuring that set containers/components do not come in contact with floor. Break both cannulae on amotosalen container.
7. Allow plasma and amotosalen solution to flow into illumination container marked by number "1".
8. Ensure that plasma drains completely from initial plasma container into illumination container.
9. Express air from the illumination container into the amotosalen container.
10. When air is removed and plasma has fully drained back into the illumination container, mix illumination container thoroughly by gentle agitation to ensure complete mixing of amotosalen and plasma.
11. Express a small amount of plasma and amotosalen mixture into tubing, filling at least 4cm of tubing. .
12. Seal tubing between illumination container and amotosalen container within the 4cm.

Warning : During illumination, tubing must be held within large compartment of illumination tray.
13. Remove and discard initial plasma container, amotosalen container and excess tubing.

C- Illumination

Illuminate plasma

Refer to INTERCEPT Illuminator Operator's Manual for complete illumination instructions for use. **Warning :** All plasma, both in illumination container and tubing, must be within large compartment of illuminator tray in order for inactivation to occur. The process requires unimpeded light transmission through tray and illumination container with plasma. No labels or other material should be on this area. Tray must be clean. Labels should be placed on illumination container flap only. Illumination container should lay flat.

D- Processing with Compound Adsorption Device (CAD)

1. Hang illumination container, allowing CAD to hang freely, with storage containers kept in an inverted position.
2. Close clamp on bypass line; ensure all other clamps are open.
3. Break cannula on illumination container and allow plasma to flow through CAD marked with the number "2" into storage containers marked with the number "3".
4. Once plasma has emptied from illumination container and passed through CAD, close clamp on tubing leading from the CAD and open clamp on the bypass line.
5. Hang CAD together with illumination container.
6. Remove storage containers from tab on CAD and allow them to hang ports up.
7. Completely express air from storage containers through bypass line.
8. Close clamp on the bypass line and open the clamps on storage containers.
9. Open clamp below the CAD, allowing plasma to drain into storage containers.
10. Ensure that the storage container tubing contains plasma and no air. Close clamps on storage containers.
11. Re-distribute plasma volume between storage containers, if necessary.
12. Ensure appropriate donor identification is attached to each filled plasma storage container.
13. Disconnect storage containers from set by heat-sealing, allowing sufficient tubing length for segments.
14. Discard CAD and illumination container. The INTERCEPT Plasma process is now complete.
15. Seal tubing as appropriate for making segments as needed.
16. Follow internal procedures for freezing plasma.

Caution: DEHP is known to be released from polyvinyl chloride (PVC) medical devices; increased leaching can occur with extended storage or increased surface area contact. The INTERCEPT processing sets only have tubing components, container ports and, if included, an in-line filter that contain PVC; all containers and other parts are PVC-free. During use of this processing set, blood components are in contact with PVC for a brief period of time (approx. <15 minutes). Based on limited surface area contact and minimal contact time, DEHP levels in blood components after use of the processing set are estimated to be well below those resulting from other medical applications containing PVC tubing (e.g. hemodialysis, intravenous fluid administration, extracorporeal membrane oxygenation and cardiopulmonary bypass procedures). The risks associated with DEHP released to the blood components must be weighed against the benefits of therapeutic transfusion and inactivation of harmful viruses, bacteria and other pathogens.

Cerus, INTERCEPT, and INTERCEPT Blood System are trademarks of Cerus Corporation

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Cerus, INTERCEPT, et INTERCEPT Blood System sont des marques de Cerus Corporation.