
**I-TROMBOSİT HAVUZLAMA SETİNİN TANIMI****Ahtapot konfigürasyonu**

REF 20304046555

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
**DEFINITION OF THE I-PLATELET POOLING SET**

**Octopus configuration**

**REF** 20304046555

One 700 ml pooling bag with 8 tubing leads for connection of buffy coats (BC) and one additional tubing lead for Platelet Additive Solution (PAS) dilution, an in-line Sepacell™ PLX-5 platelet leukodepletion filter connected to a 1,3 liter platelet storage bag with a 100 mL sampling pouch (Figure 1).


**Figure-1**

**Train configuration**

**REF** 20304146520

A tubing lead with in line Sepacell™ PLX-5 platelet leukodepletion filter connected to a 1,3 liter platelet storage container with a 100 mL sampling pouch (Figure 2).

**Figure-2**

**INDICATIONS AND USAGE**

The I-Platelet Pooling Set is intended for the preparation of leukocyte depleted platelets obtained from whole blood derived pooled BC in additive solution. Validations were conducted using BC platelets prepared from pools of up to 8 BC.

- Up to 2 therapeutic doses of leukocyte depleted platelets can be obtained by pooling 7 or 8 BC and further treatment by a Pathogen Inactivation System (INTERCEPT™ Blood System)
- One therapeutic dose of leukocyte depleted platelets can be obtained by pooling 4 to 6 BC and further treatment by a Pathogen Inactivation System (INTERCEPT™ Blood System)
- One therapeutic dose of leukocyte depleted platelets can be obtained by pooling 4 to 6 BC without further Pathogen Inactivation treatment. Platelets can be stored at 22+ 2 °C under agitation up to 24 hours in the platelet storage container of the pooling set. The platelet concentrate must be transferred to another bag with infusion ports for transfusion.

**PRECAUTIONS**

- This product is intended to be used by certified Blood Bank personnel in the Blood Bank.
- A strict standardized and validated procedure for the preparation of venipuncture site should be in operation including skin disinfection and blood diversion of the first aliquots.
- Follow validated blood component preparation procedures per local requirements.
- In order to protect patients from an increased exposure (pools of > 6 BC) it is highly recommended that the leukocyte depleted platelets obtained with the pooling set are further processed immediately with a Pathogen Inactivation System (INTERCEPT™ Blood System) before the transfusion.
- In case of further processing with a Pathogen Inactivation System, strictly follow the instructions for use of the INTERCEPT™ Blood system.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them and, they should be loaded in such a way to prevent them from falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** risk of infection when used more than once.

- For the explanation of symbols, see Symbols Description Table.

**CAUTION:** This medical device contains di(2-ethylhexyl) phthalate (DEHP). Some subpopulations of patients may be at increased risk due to exposure to DEHP following transfusions, these patients are infants undergoing exchange transfusion, children undergoing certain medical procedures and adults undergoing ECMO (extracorporeal membrane oxygenation). It has been indicated by animal experiments that nonhuman primates are less sensitive than rodents to DEHP exposure following oral administration also. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

## **PROCESS STEPS**

### **A. GENERAL PREPARATION BEFORE POOLING**

1. For optimizing the platelet yield in the platelet concentrate, pre-selection of BC is recommended (donor platelet pre-count of at least  $250 \times 10^9/L$ ).
2. Read manufacturer instruction for use for blood bags, centrifuge, blood component separator and procedures meeting local validation requirements to produce BC.
3. The recommended separation time for optimal filtration of the platelet concentrate is approximately 6 minutes. The blood component separator should be set up properly before use.

### **B. POOLING OF BUFFY COATS**

#### **Preparation before pooling**

1. Select BC previously separated and stored.
2. Open the primary package and check the I- Platelet Pooling Set for any visible defects. Do not use defective sets.

#### **Train method of pooling (indicative)**

1. Sterile connect BC one to each other using the top and bottom tubings of the bags of the train.
2. Sterile connect the first BC of the train to a Platelet Additive Solution (PAS) bag.
3. Sterile connect the last (bottom) BC bag to the collection tube of the pooling set leaving 15-20 cm of tubing and close the clamp.
4. Place a re-usable clamp on the tube between the PAS bag and the first BC bag.
5. Hang the system from the PAS bag.
6. Let all BC flow by gravity into the last bag of the train.
7. To maximize the BC content recovery, rinse each BC bag sequentially using a portion of the PAS fluid. Let the remaining volume of PAS flow into the last bag of the train. Seal and detach the rest of the train from the last bag containing the pool.

**Warning:** Make sure that outlet ports are free from red blood cells (RBC) before centrifugation. RBC may pass into the platelet storage container or block the filter.

#### **Octopus method of pooling (indicative)**

1. Sterile connect BC (up to 8) to each dedicated tubing lead of the Octopus.
2. Close all the clamps.
3. Sterile connect a PAS bag to the last tubing lead of the Octopus.
4. Hang the set.
5. Open the clamps and let all the BC flow into the pooling bag.
6. Close the clamps.
7. To maximize the BC content recovery, rinse each BC bag sequentially using a portion of the PAS fluid. Let the remaining volume of PAS flow into the pooling bag. Seal and detach the Octopus from the pooling bag.

**Warning:** The pooling bag should be filled with a maximum of 650 ml of BC and PAS. Make sure that, the outlet ports are free from RBC before centrifugation. RBC may pass into the platelet storage container or block the filter.

### **C. CENTRIFUGATION**

1. Pack the assembly of pooled buffy coats into centrifuge buckets according to blood center validated procedure.
2. Choose the appropriate platelet separation program. Follow user instructions of the centrifuge manufacturer.

**Warning:** Pack the sets into the centrifuge buckets carefully in order to avoid damages to the bags or poor separation results.

### **D. SEPARATION**

1. After centrifugation, place the set on the blood component separator.
2. Hang the pooling bag onto the pins of the separator device.
3. Ensure the leukodepletion filter, tubings and platelet storage bag are placed on the device properly. Close the slide clamp of the sampling bag.
4. Break the breakaway cannula of the pooling bag. Start the device using the appropriate platelet separation program. Follow user instructions of the manufacturer.
5. The Press plate will press the collection bag and platelets will flow into the storage bag.
6. When RBC reach the tube downstream of the filter, stop the separation.
7. Seal and detach the platelet storage bag from the rest of the assembly.
8. To remove the air, gently fold the storage bag and squeeze the air into the sampling pouch. Once the air is removed, close the clamp and seal off the storage bag allowing to keep the 10-15 cm tubing for further processing with the pathogen inactivation system.

### **E. STORAGE**

1. If further processing of the platelet concentrate is conducted with the INTERCEPT™ Blood System, platelets must be stored at  $22^\circ\text{C} \pm 2^\circ\text{C}$  under continuous agitation. Strictly follow the instructions for use of the INTERCEPT™ Blood System for further processing
2. If no further processing with a Pathogen Inactivation System will be performed, the leukocyte depleted platelets obtained with the pooling set should be used within 24 hours from the platelet pool production.

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