

en - Do not use if the product, its sterile barrier system or its packaging is damaged or shows any sign of deterioration.

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en - Sterile fluid path. Sterilized by a combination of steam and radiation.

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en - Refer to Instructions for Use

fr - Consulter la notice

en - Protect from direct sunlight and strong UVA light source.

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en - Contains or presence of phthalates: Di-(2-ethylhexyl)phthalate (DEHP)

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## INTERCEPT PROCESSING SET for SMALL VOLUME PLATELET UNITS

### For use with INTERCEPT Illuminator

Each set is wrapped in a tamper-evident package and includes one 15mL 3mM amotosalen hydrochloride solution container (Formula: Amotosalen HCl 101mg - Natr. chlorid. 924mg - Aqua ad iniecd. ad 100mL), one illumination container, one container with Compound Adsorption Device (CAD), one INTERCEPT platelet storage container. The set is sterilized by a combination of steam and radiation.

### Indications and Usage

This set is used with an INTERCEPT Illuminator to inactivate a broad spectrum of viruses, bacteria and leukocytes in platelets. INTERCEPT platelets are indicated for support of patients requiring platelet transfusions, according to clinical practice guidelines. INTERCEPT Platelets may be stored up to 7 days from time of collection, at 20-24°C with continuous agitation. INTERCEPT Platelets stored up to 7 days have been shown to adequately prevent and control bleeding. Any extension of platelet storage time from current blood center limits should be evaluated per Directive [2004/33/EC] and validated according to local blood bank procedures.

Platelet additive solutions approved for use with INTERCEPT: InterSol, SSP+

### Contraindications

Use of INTERCEPT platelets is contraindicated in patients with a history of allergic response to amotosalen or psoralens.

### Precautions

Do not use if : tamper-evident package has been opened; signs of deterioration are visible; fluid path closures are loose or not intact; cannulae are broken or there is no fluid in amotosalen solution container.

Do not store above 25°C. Do not vent. Do not freeze. Protect the pack and tubing from sharp objects.

Unused sets in open aluminium foil may be kept 20 days at room temperature by folding and securing open end of aluminium foil. Units removed from the aluminium foil must be used within 8 hours.

Keep set in light-protective package until time of use. Protect from direct sunlight and strong UVA light source.

Set is single use only. Do not reuse. Do not sterilize. This product is not designed for reuse. Misuse can result in adverse reactions, including severe illness and possibly death.

All the following conditions must be met for pathogen inactivation:

- Platelets must be prepared with additive solution for platelet re-suspension in the volume range specified in **Table 1**.
- Platelet count, volume and red blood cell (RBC) content must be within ranges specified in **Table 1**.
- Platelets mixed with amotosalen must be exposed to UVA light dose from INTERCEPT Illuminator. No other source of UVA light may be used.
- Platelets collected on Day 0 must be exposed to UVA light by end of Day 1.
- After illumination, platelets must be agitated in CAD container in accordance with the duration specified in **Table 2**.

This process is designed to be a closed system. Treatment with INTERCEPT Blood System does not replace applicable standards for processing in open and closed systems. If there is a leak in the set during processing, platelet product must be discarded.

**Warnings :** Amotosalen in contact with skin may result in photosensitization in the presence of ultraviolet light. If skin exposure occurs, flush exposed skin copiously with water. Sterile connecting device (SCD) should be used according to manufacturer's instructions for use.

### Notes to Physicians

While laboratory studies of amotosalen processing with UVA light have shown a reduction in levels of certain viruses and bacteria, there is no pathogen inactivation process that has been shown to eliminate all pathogens.

INTERCEPT platelet components should not be prescribed to newborns treated with phototherapy devices that emit a peak energy wavelength less than 425 nm, and/or whose emission bandwidth limit is less than 375 nm, due to the risk of erythema resulting from potential interaction between ultraviolet light (below 400 nm) and residual amotosalen.

### Instructions for Use

**Materials Needed :** One (1) INTERCEPT Processing Set for Small Volume Platelet Units.

**Equipment Needed :** INTERCEPT Illuminator, Sterile Connecting Device (SCD), Tube Sealer, Flatbed Agitator.

Temperatures should be controlled to meet applicable regulations for platelet processing.

Used and unused INTERCEPT sets should be discarded like any used blood containers, as biohazardous waste.

### Process Steps

#### A- Preparation of Platelets

Platelets should be prepared with additive solution. Platelets are suspended in approximately 35% plasma and 65% additive solution (**Table 1**). Platelet products within the following ranges have been shown to be acceptable for use with this processing set.

**Table 1**

Platelet Count	Volume	Plasma Content	Additive Solution Content	RBC Content
2.5 - 6.0 x 10 <sup>11</sup>	255-325mL	32-47%	53-68%	<4 x 10 <sup>6</sup> RBC/mL

#### B- Amotosalen Addition to Platelets

1. Remove set from package. Unwrap only illumination container from organizer.
2. Weld tubing from platelet container to amotosalen container tubing using SCD.
3. Label set containers using appropriate donation identification. See warning in Section C. 4. Hang platelets and break both cannulae on amotosalen container.
5. Allow platelets and amotosalen solution to flow into illumination container, marked by number "1".
6. Ensure that platelets are completely transferred to illumination container by expressing air from platelet and amotosalen container into illumination container.
7. When platelets are in illumination container, mix thoroughly by gentle agitation to ensure complete mixing of amotosalen and platelets.
8. Express air from platelets into amotosalen container.
9. Express a small amount of platelet and amotosalen mixture into tubing, filling about 4cm of tubing.
10. Seal tubing between illumination container and amotosalen container, so that tubing is no longer than approximately 4cm from illumination container inlet port.  
**Warning :** During illumination, tubing must be held within large compartment of illumination tray.
11. Remove and discard empty platelet and amotosalen containers.

#### C- Illumination

Illuminate platelets.

Refer to INTERCEPT Illuminator Operator's Manual for complete illumination instructions for use. **Warning :** Platelets in illumination container must be within large compartment of illuminator tray in order for inactivation to occur. The process requires unimpeded light transmission through tray and illumination container with platelets. No labels or other material should be on this area. Tray must be clean. Labels should be placed on illumination container flap only. Illumination container should lay flat.

#### D- Agitation with CAD

**Warning :** Do not fold or bend CAD.

1. Unwrap set from organizer.
2. Hang platelets, break cannula and allow platelets to flow into CAD container, marked by number "2".
3. Express air from CAD container into illumination container.
4. Seal tubing close to inlet port of CAD container.
5. Remove and discard empty illumination container.
6. Place CAD container on flatbed agitator for the duration specified in **Table 2**.

**Table 2**

Suspension Medium	CAD Agitation Duration
Additive Solution (53 - 68%)	4 to 16 hours

#### E- Transfer to Storage Container

1. Remove platelets from agitator and hang platelets.
2. Close clamp on platelet sampling pouch.
3. Break cannula and allow platelets to flow into storage container, marked by number "3".
4. Express air from storage container into CAD container.
5. Seal tubing close to inlet port of storage container.
6. Remove and discard empty CAD container. Place platelets on platelet agitator.

The INTERCEPT Platelet process is now complete.

#### F- Sampling Platelet Product (optional)

1. Mix INTERCEPT platelets by gently agitating platelet storage container.
2. Open clamp to platelet sampling pouch and squeeze several times.
3. Allow sampling pouch to fill with platelets. Seal tubing.
4. Remove sample pouch.
5. Transfer sample to appropriate laboratory tube immediately.

**Caution:** DEHP is known to be released from polyvinyl chloride (PVC) medical devices; increased leaching can occur with extended storage or increased surface area contact. The INTERCEPT processing sets only have tubing components, container ports and, if included, an in-line filter that contain PVC; all containers and other parts are PVC-free. During use of this processing set, blood components are in contact with PVC for a brief period of time (approx. <15 minutes). Based on limited surface area contact and minimal contact time, DEHP levels in blood components after use of the processing set are estimated to be well below those resulting from other medical applications containing PVC tubing (e.g. hemodialysis, intravenous fluid administration, extracorporeal membrane oxygenation and cardiopulmonary bypass procedures). The risks associated with DEHP released to the blood components must be weighed against the benefits of therapeutic transfusion and inactivation of harmful viruses, bacteria and other pathogens.

Cerus, INTERCEPT, and INTERCEPT Blood System are trademarks of Cerus Corporation  
InterSol is a trademark of Fenwal Inc.  
SSP+ is a trademark of MacoPharma

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Cerus, INTERCEPT, et INTERCEPT Blood System sont des marques de Cerus Corporation.
InterSol est une marque de Fenwal, Inc. SSP+ est une marque de MacoPharma